

# ome\_demo

June 8, 2021

## 1 ezomero Demo

## 2 OME Community Meeting 2021

### 2.1 Before we begin...

For this demo, I am connecting to a JAX OMERO server via our VPN, so you will not be able to follow along exactly.

I am using a MacBook Pro running macOS Catalina, working in a simple conda environment in which I installed:

- python 3.7
- ezomero (via pip)
- jupyterlab (via conda-forge)
- matplotlib

All other dependencies (including OMERO.py) were automatically installed.

### 2.2 Create a connection

#### 2.2.1 Create a connection from parameters

```
[1]: # All functions in ezomero have useful docstrings
import ezomero

%pdoc ezomero.connect
```

**Class docstring:**

Create an OMERO connection

This function will create an OMERO connection by populating certain parameters for ``omero.gateway.BlitzGateway`` initialization by the procedure described in the notes below. Note that this function may ask for user input, so be cautious if using in the context of a script.

Finally, don't forget to close the connection ``conn.close()`` when it is no longer needed!

Parameters

```

-----
user : str, optional
      OMERO username.

password : str, optional
          OMERO password.

group : str, optional
       OMERO group.

host : str, optional
      OMERO.server host.

port : int, optional
      OMERO port.

secure : boolean, optional
        Whether to create a secure session.

config_path : str, optional
             Path to directory containing '.ezomero' file that stores connection
             information. If left as ``None``, defaults to the home directory as
             determined by Python's ``pathlib``.

```

#### Returns

```

-----
conn : ``omero.gateway.BlitzGateway`` object or None
      OMERO connection, if successful. Otherwise an error is logged and
      returns None.

```

#### Notes

```

-----
The procedure for choosing parameters for ``omero.gateway.BlitzGateway``
initialization is as follows:

```

- 1) Any parameters given to `ezconnect` will be used to initialize  
`omero.gateway.BlitzGateway`
- 2) If a parameter is not given to `ezconnect`, populate from variables  
in ``os.environ``:
  - \* OMERO\_USER
  - \* OMERO\_PASS
  - \* OMERO\_GROUP
  - \* OMERO\_HOST
  - \* OMERO\_PORT
  - \* OMERO\_SECURE
- 3) If environment variables are not set, try to load from a config file.

This file should be called '.ezomero'. By default, this function will look in the home directory, but ``config\_path`` can be used to specify a directory in which to look for '.ezomero'.

The function ``ezomero.store\_connection\_params`` can be used to create the '.ezomero' file.

Note that passwords can not be loaded from the '.ezomero' file. This is to discourage storing credentials in a file as cleartext.

- 4) If any remaining parameters have not been set by the above steps, the user is prompted to enter a value for each unset parameter.

Call docstring:

Call self as a function.

```
[2]: conn = ezomero.connect(user='djme',
                           group='Research IT',
                           host='bhomero01lp',
                           port=4064,
                           secure=True)
```

Enter password: .....

```
[3]: ezomero.print_groups(conn)
```

Groups:

```
          default: 3    member
          Public: 53
Microscopy_Service: 103
          Canine_PDX: 104
          Murray Lab: 105
          Korstanjelab: 106
          MTB: 107
          KOMP_eye: 153  member
          Hinson Lab: 203
          Nishina Lab: 204
          JAX_Hacks: 253  owner
Bolcun-Filas Lab: 303
          Research IT: 304  owner
          Verhaak Lab: 353  owner
          KOMP_histopath: 403  member
          Cube: 453
          OConnell Lab: 503
          Robson Lab: 504
          Braun Lab: 505
          Rosenthal Lab: 553
          GRS: 603
```

```
[4]: print(conn.group)
```

Research IT

```
[5]: conn.close()
```

### 2.2.2 Create a connection from prompts

```
[6]: conn = ezomero.connect()
```

```
Enter username: djme
Enter password: .....
Enter group name (or leave blank for default group):
Enter host: bhomero01lp
Enter port: 4064
Secure session (True or False): t
```

```
[7]: conn.close()
```

### 2.2.3 Create a connection from stored parameters

```
[8]: ezomero.store_connection_params()
```

```
Enter username: djme
Enter group name (or leave blank for default group):
Enter host: bhomero01lp
Enter port: 4064
Secure session (True or False): t

Connection settings saved to /Users/djme/.ezomero
```

```
[9]: # Note you still get prompted for password
     conn = ezomero.connect()
```

```
Enter password: .....
```

```
[10]: conn.close()
```

### 2.2.4 Create a connection from environment variables

```
[11]: # first remove the stored connection parameters

import os
try:
    os.remove('/Users/djme/.ezomero')
except FileNotFoundError:
    pass
```

```
# Setting these using some jupyter magic.
# You wouldn't necessarily do this in practice
%env OMERO_USER=djme
%env OMERO_GROUP=JAX_Hacks
%env OMERO_HOST=bhomer01lp
```

```
env: OMERO_USER=djme
env: OMERO_GROUP=JAX_Hacks
env: OMERO_HOST=bhomer01lp
```

```
[12]: # Since we didn't use all environment variables, we get prompts
      conn = ezomero.connect()
```

```
Enter password: .....
Enter port: 4064
Secure session (True or False): True
```

```
[13]: conn.group
```

```
[13]: 'JAX_Hacks'
```

```
[14]: conn.close()
```

## 2.3 Browse OMERO data

```
[15]: conn = ezomero.connect('djme',
                             group='Research IT',
                             host='bhomer01lp',
                             port=4064,
                             secure=True)
```

```
Enter password: .....
```

```
[16]: ezomero.print_projects(conn)
```

```
Projects:
  test_friday: 501
  zarr: 553
  testdata: 559
  test: 652
  test_auto: 902
  test_auto2: 903
  test_auto1: 904
```

```
[17]: ezomero.print_datasets(conn, 501)
```

```
Datasets in Project "test_friday":
  test_ds_friday: 1456
```

```
[18]: im_ids = ezomero.get_image_ids(conn, dataset=1456)
      print(im_ids)
```

```
[24902, 24901]
```

```
[19]: import numpy as np


      # grab an image
      im_obj, pix = ezomero.get_image(conn, 24901)
      pix.shape # tzyxc (skimage convention)
```

```
[19]: (1, 1, 768, 768, 3)
```

```
[20]: im = np.squeeze(pix)
      im.shape
```

```
[20]: (768, 768, 3)
```

```
[21]: import matplotlib.pyplot as plt
      plt.imshow(im);
```



## 2.4 Creating Projects, Datasets, and Images

```
[22]: ezomero.print_projects(conn)
```

```
Projects:
    test_friday:    501
    zarr:    553
    testdata:    559
    test:    652
    test_auto:    902
    test_auto2:    903
    test_auto1:    904
```

```
[23]: pid = ezomero.post_project(conn, "OMEdemo2021", description='This is a demo')
      print(pid)
```

913

```
[24]: ezomero.print_projects(conn)
```

```
Projects:
    test_friday:    501
    zarr:    553
    testdata:    559
    test:    652
    test_auto:    902
    test_auto2:    903
    test_auto1:    904
    OMEdemo2021:    913
```


```
[25]: did = ezomero.post_dataset(conn, "test dataset", description='This is our_
      ↪dataset', project_id=pid)
      print(did)
```

2586

```
[26]: ezomero.print_datasets(conn, pid)
```

```
Datasets in Project "OMEdemo2021":
    test dataset:    2586
```

```
[27]: # Crop image
      crop_im = im[:300,:300,:]
      plt.imshow(crop_im);
```



```
[28]: crop_im.shape
```

```
[28]: (300, 300, 3)
```

```
[29]: # Need to rearrange axes into OMERO style (XYZCT)
```

```
omero_im = np.moveaxis(crop_im[:, :, None, :, None], 0, 1)
omero_im.shape
```

```
[29]: (300, 300, 1, 3, 1)
```

```
[30]: im_id = ezomero.post_image(conn, omero_im,
                                image_name='gradient_image',
                                dataset_id=did,
                                source_image_id=24901)
ezomero.get_image_ids(conn, did)
```

WARNING:root:Using this function to save images to OMERO is not recommended when `transfer=ln\_s` is the primary mechanism for data import on your OMERO instance. Please consult with your OMERO administrator.

```
[30]: [128137]
```



## 2.5 Working with MapAnnotations

```
[31]: # The key-value pairs of MapAnnotations are basically dicts
      d = {'key1': 'value1',
           'antibody': 'my favorite ab',
           'temp (C)': 27}
      ns = 'jax.org/example/namespace'

      map_ann_id = ezomero.post_map_annotation(conn, 'Image', im_id, kv_dict=d, ns=ns)
      print(map_ann_id)
```

102287

```
[32]: # Retrieve the MapAnnotation we just posted
      map_ann = ezomero.get_map_annotation(conn, map_ann_id)
      print(map_ann)
```

{'key1': 'value1', 'antibody': 'my favorite ab', 'temp (C)': '27'}

```
[33]: # Change a value and update original MapAnnotation
      map_ann['temp (C)'] = 25
      ezomero.put_map_annotation(conn, map_ann_id, map_ann)
      new_map_ann = ezomero.get_map_annotation(conn, map_ann_id)
      print(new_map_ann)
```

{'key1': 'value1', 'antibody': 'my favorite ab', 'temp (C)': '25'}

```
[34]: conn.close()
```